



Guardian Wireless Alarm System



Product Code: HJM340

Description:

The Guardian is a modern wireless falls alarm system. Its fundamental purpose is to act as a pressure-sensitive warning system to alert nurses and carers to unexpected and potentially dangerous movements from users. Triggered by a release or applied pressure to a surface, the Guardian can be configured to suit an extensive range of situations and requirements.

Key Features:

- Safe and versatile wireless system
- Compatible with most nurse call system*
- Simple to set up and use
- Up to 12 months of battery life for the transmitter and up to 3 months for the receiver.
- Pair with up to 8 transmitters.
- Can be paired with a variety of sensor mats to suit different situations.
- Automatically adjusts the trigger depending on which sensor mat is used.
- Patient care function reverts to normal operation after 10 minutes (no alarm triggers while giving care). Works on both normally open and normally closed systems.
- Silent alarm function (activates nurse call system with no noise at the source).
- Low battery warning for receiver and transmitters (before function is interrupted).

Additional Products:

- Power Supply [HJR079](#)
- Seat Mat with Pocket [HJM343S](#)
- Mattress Mat with Pocket [HJM343M](#)
- Floor Mat [HJM301](#)
- Adapter Cable to use HLV301S and HLV301M [HJM342](#) Transmitter [HJM341](#)

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Guardian is guaranteed for 1 year from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk;

w: www.harvesthealthcare.co.uk

Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: August 2023